



Service is our passion.

WKD or accredited calibration?

WKD: Factory calibration of environmental simulation systems

- Factory calibrations may be performed by any calibration laboratory, no accreditation / certification is required for this.
- WKD is a method developed by weisstechnik[®] itself procedure.
- It can be carried out with any measuring equipment that is calibrated and meteorologically traceable.
- The layout of the calibration certificate is determined by the calibration laboratory itself.
- WKD provides accurate results, but is not an internationally accepted procedure.



Accredited calibration of environmental simulation systems

- Execution only by accredited bodies that calibrate according to calibrate according to DKD-R 5-7.
- The weisstechnik[®] calibration laboratory is accredited by the DAkkS according to DIN EN ISO/IEC 17025.
- The used measuring equipment must be traceable to national standards.
- The layout of the calibration certificate is firmly defined and must be adhered to.
- The measurement results are internationally comparable.
- An accredited calibration offers the highest possible reliability.



We are there for you 365 days a year, 24 hours a day.

Our service technicians are excellently trained specialists. Our service teams are professionally competent and available at all times to provide you with reliable support in a spirit of partnership.

More safety from a single source.

Regardless of which calibration option you choose, we can adjust your system promptly in the event of deviations. You can also avoid unnecessary downtime by calibrating your system as part of annual maintenance.

In future, take advantage of our fast, uncomplicated and, above all, more cost-effective way of having your systems serviced, calibrated and adjusted.

We ensure trouble-free, calibrated and legally compliant operation of your systems!

Request a non-binding offer now:
service-vertrieb.ut@weiss-technik.com



We measure ourselves by service!

**Our services -
plenty of good arguments:**

- Global service network
- Wide range of preventive maintenance services
- Reliable spare parts supply
- Special assignments available at any time
- Professional disposal of your old appliances with proof

Our service experts are always close by.

**24/7-Service-Helpline:
+49 1805 666 556**

Weiss Technik GmbH
Greizer Straße 41-49
35447 Reiskirchen/Germany
T +49 6408 84-0
info@weiss-technik.com
weiss-technik.com



Test it. Heat it. Cool it.