

Laboratory Oven

Laboratory oven	UT 6060	T 6030	UT 6420
Test space volume (l)	52	30	375
Test space dimensions (WxDxH)	405 x 340 x 380mm	352 x 370 x 231mm	544 x 522 x 1319mm
External dimensions (WxDxH)	745 x 535 x 555mm	552 x 535 x 552mm	744 x 715 x 1707mm
Temperature range	up to 300°C	up to 300°C	up to 300°C
Humidity range	-		
Dew point range	/		
Type cooling	/		
Tension	1 x 230 V + PE	1 x 230V + PE	3 x 400V + N + PE
Current	16 A	3,4 A	16 A

Laboratory oven	VT 6060 M
Serial number	93111268
Year of construction	1993
Test space volume (l)	53
Test space dimensions (WxDxH)	415 x 345 x 371mm
External dimensions (WxDxH)	744 x 630 x 576mm
Temperature range	Up to 200°C
Tension	1 x 230 V + PE
Current	7,4 A

Laboratory oven	OMS 60
Test space volume (l)	66l
Test space dimensions (WxDxH)	354x368x508
External dimensions (WxDxH)	530x565x720 - 42kg
Temperature range	+10 K tot +250°C
Tension	230V 1/PE AC 50/60 Hz

Laboratory oven	OMH 750-3P
Test space volume (l)	731l
Test space dimensions (WxDxH)	1004x545x1335
External dimensions (WxDxH)	1261x770x1545 - 185kg
Temperature range	+10 K tot +300°C
Tension	400V 3/PN AC 50/60 Az